

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.provita.nt-rt.ru | pia@nt-rt.ru

Технические характеристики на медицинские стационарные сборные системы Direct, Direct Plus компании **PROVITA**

Виды товаров: держатели для защиты штекера, протекторы настенного крепления, секции, настенные блоки с трубками, панели, протекторы для гигиенической задней панели, шкафные устройства для крепления ультразвуковых датчиков, базовые держатели для горизонтальных и вертикальных настенных направляющих, базовые держатели для круглой трубы, держатели чреспищеводного зонда для операционных столов и др.



Product highlights

WALL MOUNT PROTECTOR, GREEN INSERT 1.200 mm


Product number: T-W72416

working length:	1200 mm
Field of application:	TEE-Care (Ultrasound probes)
Function of the unit:	storage of one probe
Material of housing:	Aluminium, pulverbeschichtet
Secure Support:	Cable deflection roller for optimal cable routing
Marking insert:	green = endoscope reprocessed
Swing guard:	up to Ø 23 mm
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Total load capacity (kg):	2

Packaging dimensions

Packaging unit: Unit

Width Box1	0,15 m
Height Box1	0,1 m
Length Box1	1,75 m
Weight Box1	5,5 kg



Product highlights

WALL MOUNT PROTECTOR, RED INSERT 1.200 mm

Product number:	T-W72426
Field of application:	TEE-Care (Ultrasound probes)
Working length:	1.200
Function of the unit:	storage of one probe
Material of housing:	Aluminium, pulverbeschichtet
Secure Support:	Cable deflection roller for optimal cable routing
Marking insert:	red = endoscope contaminated
Swing guard:	up to Ø 23 mm
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Total load capacity (kg):	2

Packaging dimensions

Packaging unit: Unit

Width Box1	0,15 m
Height Box1	0,1 m
Length Box1	1,75 m
Weight Box1	5 kg



Product highlights

HOLDER FOR PLUG PROTECTION

Wall mount

Product number: E-B-S006

Material: Stainless steel

Size: 150 x 55 mm

Packaging dimensions


Packaging unit: Unit

Width Box1 0,1 m

Height Box1 0,1 m

Length Box1 0,2 m

Weight Box1 0,5 kg



Product highlights

Application

WALL UNIT "DIRECT+" probe attachment

Product number: T-W07401

Area of application:	TEE-Care
Tray material:	Stainless steel
Function of the unit:	storage of 1 probe
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by wipe disinfection
Size:	appr. 400 x 1.700 mm

Packaging dimensions

Packaging unit: Unit

Width Box1	0,45 m
Height Box1	0,1 m
Length Box1	1,75 m
Weight Box1	8 kg



Product highlights



TIME-SAVING CLEANING

Application

WALL UNIT "DIRECT+" WITH 1 TUBE immersion disinfection

Product number: T-T07432

Cleaning/disinfection tubes :	Immersion disinfection
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Function of the unit:	desinfection without pre-cleaning
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by immersion disinfection
Type of Fixation:	To screw on
Size:	ca. 240 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,29 m
Height Box1	0,15 m
Length Box1	1,75 m
Weight Box1	6,5 kg



Product highlights

WALL UNIT "DIRECT+" WITH 1 TUBE leak test

Product number: T-M07432

Cleaning/disinfection tubes :	leak test
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Function of the unit:	storage of 1 probe
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Size:	ca. 240 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,29 m
Height Box1	0,15 m
Length Box1	1,75 m
Weight Box1	6 kg



Product highlights

Application

WALL UNIT "DIRECT+" WITH 1 TUBE

Immersion disinfection, optional protector


Product number: T-T08436

Cleaning/disinfection tubes :	Immersion disinfection
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Disinfection:	Immersion disinfection
Function of the unit:	desinfection without pre-cleaning
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by immersion disinfection
Type of Fixation:	To screw on
Size:	ca. 360 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,41 m
Height Box1	0,15 m
Length Box1	1,75 m
Weight Box1	7,5 kg



Product highlights

WALL UNIT "DIRECT+" WITH 2 TUBES pre-cleaning and immersion disinfection

Product number: T-T07466

Cleaning/disinfection tubes :	Cleaning and Immersion disinfection
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Disinfection:	Immersion disinfection
Function of the unit:	disinfection including pre-cleaning
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by immersion disinfection
Type of Fixation:	To screw on
Size:	ca. 360 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,41 m
Height Box1	0,15 m
Length Box1	1,75 m
Weight Box1	10 kg



Product highlights

WALL UNIT "DIRECT+" WITH 2 TUBES immersion disinfection and leak test


Product number: T-T07456

Cleaning/disinfection tubes :	Immersion disinfection and Leak test
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Function of the unit:	desinfection without pre-cleaning
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by immersion disinfection
Type of Fixation:	To screw on
Size:	ca. 360 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,41 m
Height Box1	0,15 m
Length Box1	1,75 m
Weight Box1	8,5 kg



Product highlights

WALL UNIT "DIRECT+" WITH 2 TUBES protector optional


Product number: T-T08466

Cleaning/disinfection tubes :	Cleaning and Immersion disinfection
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Function of the unit:	disinfection including pre-cleaning
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by immersion disinfection
Type of Fixation:	To screw on
Size:	ca. 445 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,45 m
Height Box1	0,15 m
Length Box1	1,75 m
Weight Box1	10 kg



Product highlights

WALL UNIT "DIRECT+" WITH 3 TUBES pre-cleaning, immersion des., leak test

Product number: T-T07476

Cleaning/disinfection tubes :	Cleaning, Immersion disinfection and Leak test
Field of application:	TEE-Care (Ultrasound probes)
Dip tube length:	1.160 (ca.)
Tray material:	Stainless steel
Function of the unit:	disinfection including pre-cleaning
Secure Support:	Cable deflection roller for optimal cable routing
Reprocessing of the probe:	by immersion disinfection
Type of Fixation:	To screw on
Size:	ca. 445 x 1.700

Packaging dimensions

Packaging unit: Unit

Width Box1	0,5 m
Height Box1	0,2 m
Length Box1	1,75 m
Weight Box1	11 kg



Product highlights



USER-FRIENDLY CLEANING

PROTECTOR FOR STATIONARY "DIRECT +", GREEN INSERT 1.200 mm

Product number: T-W73416

Premounted :	Hygienic rear panel
Working length:	1.200
Disinfection:	Wipe disinfection
Function of the unit:	storage of one probe
Material of housing:	Aluminium, pulverbeschichtet
Marking insert:	green = endoscope reprocessed
Swing guard:	up to Ø 23 mm
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Total load capacity (kg):	2



Product highlights

PROTECTOR FOR HYGIENIC REAR PANEL 1.200 mm

Product number: T-W73426

Premounted :	Hygienic rear panel
working length:	1200 mm
Disinfection:	Wipe disinfection
Function of the unit:	storage of one probe
Material of housing:	Aluminium, pulverbeschichtet
Marking insert:	red = endoscope contaminated
Swing guard:	up to Ø 23 mm
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Total load capacity (kg):	2

Packaging dimensions

Packaging unit: Unit

Width Box1	0,15 m
Height Box1	0,1 m
Length Box1	1,25 m
Weight Box1	5 kg



Product highlights

Application

CABINET UNIT for 1 ultrasound probe


Product number: T-W06401

Field of application:	TEE-Care (Ultrasound probes)
Tray material:	Stainless steel
Function of the unit:	storage of 1 probe
Secure Support:	Cable deflection roller for optimal cable routing
Swing guard:	up to Ø 23 mm
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Size:	ca. 200 x 600 mm

Packaging dimensions

Packaging unit: Unit

Width Box1	0,25 m
Height Box1	0,85 m
Length Box1	0,65 m
Weight Box1	2 kg



Product highlights

CABINET UNIT for 2 ultrasound probes


Product number: **T-W06402**

Field of application:	TEE-Care (Ultrasound probes)
Tray material:	Stainless steel
Function of the unit:	storage of 2 probes
Secure Support:	Cable deflection roller for optimal cable routing
Swing guard:	up to Ø 23 mm
Reprocessing of the probe:	by wipe disinfection
Type of Fixation:	To screw on
Size:	ca. 400 x 600 mm

Packaging dimensions

Packaging unit: Unit

Width Box1	0,45 m
Height Box1	0,1 m
Length Box1	0,65 m
Weight Box1	3,5 kg



Product highlights

TEE BASIC HOLDER for horizontal wall rails

Product number: T-W08010

Click system :	optional plug protection
Material:	Stainless steel/ Plastic
Plug holder:	optional against surcharge
For attachment to :	Rundrohr Ø 25 - 30 mm
Securing:	with safety bolt

Packaging dimensions

Packaging unit: Unit

Width Box1	0,3 m
Height Box1	0,2 m
Length Box1	0,3 m
Weight Box1	1,1 kg



Product highlights

TEE BASIC HOLDER for vertical wall rails and round tube

Product number: T-W08020

Click system :	optional plug protection
Material:	Stainless steel/ Plastic
Plug holder:	optional against surcharge
Securing:	with safety bolt

Packaging dimensions

Packaging unit: Unit

Width Box1	0,2 m
Height Box1	0,1 m
Length Box1	0,3 m
Weight Box1	1,2 kg



Product highlights

TEE PROBE HOLDER for operating tables

Product number: T-W08000

Click system :	optional plug protection
Probe mount:	premounted
Material:	Stainless steel/ Plastic
Plug holder:	optional against surcharge
Securing:	with safety bolt

Packaging dimensions

Packaging unit: Unit

Width Box1	0,26 m
Height Box1	0,17 m
Length Box1	0,35 m
Weight Box1	1,5 kg

Алматы (7273)495-231
Ангарск (3955)60-70-56
Архангельск (8182)63-90-72
Астрахань (8512)99-46-04
Барнаул (3852)73-04-60
Белгород (4722)40-23-64
Благовещенск (4162)22-76-07
Брянск (4832)59-03-52
Владивосток (423)249-28-31
Владикавказ (8672)28-90-48
Владимир (4922)49-43-18
Волгоград (844)278-03-48
Вологда (8172)26-41-59
Воронеж (473)204-51-73
Екатеринбург (343)384-55-89

Иваново (4932)77-34-06
Ижевск (3412)26-03-58
Иркутск (395)279-98-46
Казань (843)206-01-48
Калининград (4012)72-03-81
Калуга (4842)92-23-67
Кемерово (3842)65-04-62
Киров (8332)68-02-04
Коломна (4966)23-41-49
Кострома (4942)77-07-48
Краснодар (861)203-40-90
Красноярск (391)204-63-61
Курск (4712)77-13-04
Курган (3522)50-90-47
Липецк (4742)52-20-81

Магнитогорск (3519)55-03-13
Москва (495)268-04-70
Мурманск (8152)59-64-93
Набережные Челны (8552)20-53-41
Нижний Новгород (831)429-08-12
Новокузнецк (3843)20-46-81
Ноябрьск (3496)41-32-12
Новосибирск (383)227-86-73
Омск (3812)21-46-40
Орел (4862)44-53-42
Оренбург (3532)37-68-04
Пенза (8412)22-31-16
Петрозаводск (8142)55-98-37
Псков (8112)59-10-37

Пермь (342)205-81-47
Ростов-на-Дону (863)308-18-15
Рязань (4912)46-61-64
Самара (846)206-03-16
Саранск (8342)22-96-24
Санкт-Петербург (812)309-46-40
Саратов (845)249-38-78
Севастополь (8692)22-31-93
Симферополь (3652)67-13-56
Смоленск (4812)29-41-54
Сочи (862)225-72-31
Ставрополь (8652)20-65-13
Сургут (3462)77-98-35
Сыктывкар (8212)25-95-17
Тамбов (4752)50-40-97

Тверь (4822)63-31-35
Тольятти (8482)63-91-07
Томск (3822)98-41-53
Тула (4872)33-79-87
Тюмень (3452)66-21-18
Ульяновск (8422)24-23-59
Улан-Удэ (3012)59-97-51
Уфа (347)229-48-12
Хабаровск (4212)92-98-04
Чебоксары (8352)28-53-07
Челябинск (351)202-03-61
Череповец (8202)49-02-64
Чита (3022)38-34-83
Якутск (4112)23-90-97
Ярославль (4852)69-52-93

Россия +7(495)268-04-70

Казахстан +7(7172)727-132

Киргизия +996(312)96-26-47

www.provita.nt-rt.ru | | pia@nt-rt.ru